


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium methanesulfonate - 98%

Product Number: 304506
CAS Number: 2386-57-4
MDL: MFCD00064389
Formula: CH₃NaO₃S
Formula Weight: 118.09 g/mol



| TEST | Specification |
|---|--------------------------|
| Appearance (Color) | White |
| Appearance (Form) | Conforms to Requirements |
| Powder to Crystalline Powder to Fibers and/or Lumps and Chips | |
| Infrared spectrum | Conforms to Structure |
| Carbon | 9.8 - 10.5 % |
| Water (by Karl Fischer) | ≤ 2.5 % |

Specification: PRD.0.ZQ5.10000101676

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.