


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,3,7,8,12,13,17,18-Octaethyl-21H,23H-porphine cobalt(II)

Product Number: 257605
CAS Number: 17632-19-8
MDL: MFCD00012148
Formula: C₃₆H₄₄CoN₄
Formula Weight: 591.69 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

Purple

Appearance (Form)

Powder or Fibers

Infrared spectrum

Conforms to Structure

Wavelength

390 - 396 nm

0.002g/l, Toluene

Extinction Coefficient

≥ 135000

Specification: PRD.0.ZQ5.10000036931

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.