

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Methyl β -D-xylopyranoside - $\geq 99\%$ (GC)

Product Number:

M5878

CAS Number:

612-05-5

MDL:


MFCD00047532

Formula:

C₆H₁₂O₅

Formula Weight:

164.16 g/mol

**TEST****Specification**

Appearance (Color)

White to Off-White

Appearance (Form)

Powder

Proton NMR spectrum

Conforms to Structure

Specific Rotation

-65.0 - -60.0 °

C = 10 in H₂O at 20 deg C

Purity (GC)

 $\geq 99\%$

Solubility (Color)

Colorless to Light Yellow

Solubility (Turbidity)

Clear to Hazy

100 mg/mL, H₂O

Specification: PRD.2.ZQ5.10000010186

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.