


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Cesium formate - 98%

Product Number: 433764
CAS Number: 3495-36-1
MDL: MFCD00039103
Formula: CHCsO₂
Formula Weight: 177.92 g/mol



| TEST | Specification |
|-----------------------------|-------------------------------|
| Appearance (Color) | White |
| Appearance (Form) | Powder or Crystal or Chunk(s) |
| Infrared spectrum | Conforms to Structure |
| Purity (Titration by HCL04) | 97.5 - 102.5 % |

Specification: PRD.0.ZQ5.10000027957

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.