

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

LB Agar Ampicillin-100, Plates - pre-poured agar plates with 100 µg/ml ampicillin

Product Number:

L5667

Storage Temperature:

2 - 8 °C

TEST	Specification
------	---------------

Appearance Amber agar plates	Conforms
Sterility	Pass
Performance Testing Growth promotion and Antibiotic or Amino acid sensitivity	Pass

Specification: PRD.0.ZQ5.10000048683

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.