


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Fmoc-(S)-propargyl-Ala-OHProduct Number: **F5437**
CAS Number: 1198791-65-9Formula: C₂₁H₁₉NO₄
Formula Weight: 349.38 g/mol
Storage Temperature: -20 °C

TEST	Specification
Appearance White solid or yellow oil	Conforms
Identity H NMR or MS	Conforms
HPLC	≥ 97 %

Specification: PRD.0.ZQ5.10000063884

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.