

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com Email USA: techsery@sial.com Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal ANTI-FLAG* M2-Peroxidase (HRP) antibody produced in mouse - clone M2, purified immunoglobulin, buffered aqueous glycerol solution

Product Number:

Storage Temperature: -20 ℃

TEST Specification

Appearance (Turbidity) Clear

Colorless to Amber Appearance (Color)

Appearance (Form) Liquid Protein Conforms

by calculation: 1.1 mg/ml

ELISA Titer > 20000

Dot Blot Limit of Dection 8 ng

Minimum detection of FLAG-Fusion protein (N-BAP) when using 1:100 dilution of conjugate Dot Blot

Limit of Dection 32 ng

Minimum detection of FLAG-Fusion protein (N-BAP) when using 1:1000 dilution of conjugate Immunocytochemistry Minimum Dilution 1:100

Minimum titer of 1:100 on FLAG transfected COS cells

Specification: PRD.1.ZQ5.10000064452

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.