

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Titanium(IV) chloride solution - 0.09 M in 20% HCl

Product Number:

404985

TiCl₄

MDL:

MFCD00011267

Formula:

Cl₄Ti

Formula Weight:

189.68 g/mol

TEST	Specification
------	---------------

Appearance (Color)

Colorless

Appearance (Form)

Liquid

Mass Balance

Specification: Approximately 0.09 M TiCl₄

Purity (Titration by NaOH)

18 - 22 %

% HCl

Specification: PRD.0.ZQ5.10000027271

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.