

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:


Poly(ethylene glycol) diacrylate - average  $M_n$  1,000, contains MEHQ as inhibitor

Product Number:

729086

CAS Number:

26570-48-9



Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Color)

White to Beige

Appearance (Form)

Conforms to Requirements

Powder or Solid to Moist Solid

Proton NMR spectrum

Conforms to Structure

Molecular Number

900 - 1100

by GPC (RI, THF, PEO)

Poly Dispersity MW/MN

≤ 1.1

Water (by Karl Fischer)

Specification: PRD.0.ZQ5.10000049542

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.