

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Palladium on calcium carbonate - Evonik Noblyst® P8059 5% Pd

Product Number:

908134

**Pd**

TEST	Specification
Appearance (Color)	Grey to Dark Grey
Appearance (Form)	Powder
ICP Major Analysis	Confirmed
Confirms Palladium and Calcium Components	
Palladium (Pd)	4.5 - 5.5 % wt
Calculated to dry basis	
Activity	
Cinnamic acid hydrogenation ml/min.g	
Specification: Report Only	
Product of Supplier	Confirmed
Evonik	
Registered Trademark	Note
Noblyst is a registered trademark of Evonik Degussa GmbH	
Recommended Retest Period	-----
2 Years	

Specification: PRD.0.ZQ5.10000100558

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.