


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Water - Molecular Biology Reagent

Product Number: W4502
CAS Number: 7732-18-5
MDL: MFCD00011332
Formula: H₂O
Formula Weight: 18.02 g/mol



TEST	Specification
------	---------------

Appearance (Turbidity)	Clear
Appearance (Color)	Colorless
Appearance (Form)	Liquid
DNase, Exonuclease Detection	None Detected
NICKase, Endonuclease Detection	None Detected
RNAse Detection	None Detected
Protease Detection by FITC Casein	None Detected

Specification: PRD.1.ZQ5.10000054299

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.