

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Casein from bovine milk - technical grade

Product Number:

C7078

CAS Number:

9000-71-9

MDL:

MFCD00081481

TEST	Specification
Appearance (Color)	Off White
Appearance (Form)	Powder
% Loss on Drying	≤ 12
Nitrogen	13.5 - 15.0 %

Specification: PRD.0.ZQ5.10000037906

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.