


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:  
Nicotinic acid N-oxide - 99%

**Product Number:** 106194  
**CAS Number:** 2398-81-4  
**MDL:** MFCD00006201  
**Formula:** C<sub>6</sub>H<sub>5</sub>NO<sub>3</sub>  
**Formula Weight:** 139.11 g/mol



### TEST

### Specification

Appearance (Color)	White to Pale Yellow
Appearance (Form)	Conforms to Requirements
Powder or Crystalline Powder	
Infrared Spectrum	Conforms to Structure
Purity (Titration by NaOH)	98.5 - 101.5 %

Specification: PRD.0.ZQ5.10000040892

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.