

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Sodium hypochlorite solution - reagent grade, available chlorine 4.00-4.99 %

Product Number: 239305
CAS Number: 7681-52-9
MDL: MFCD00011120
Formula: ClNaO
Formula Weight: 74.44 g/mol
Storage Temperature: 2 - 8 °C

NaOCl

TEST	Specification
------	---------------

Appearance (Clarity)	Clear
Appearance (Color)	Colorless to Yellow and Colorless to Green-Yellow
Appearance (Form)	Liquid
Titration with Na ₂ S ₂ O ₃ Available Chlorine	4.00 - 4.99 %
Recommended Retest Period 12 Months	-----

Specification: PRD.2.ZQ5.10000016424

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.