


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium iodide dihydrate - ≥99%

Product Number: 217301
CAS Number: 13517-06-1
MDL: MFCD00149211
Formula: INa · 2H₂O
Formula Weight: 185.92 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Crystals
ICP MAJOR Analysis Confirms Na Component	Confirmed
Titration by AgNO ₃ % I	67.6 - 68.9 %

Specification Date : 04/26/2011

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.