

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Clostridium difficile Toxin A - ≥95% (SDS-PAGE)

Product Number: **SML1154**

Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Formulation	Conforms
Lyophilized from solution of 0.05 M HEPES, 0.15 M NaCl, and 5% Sucrose	
Purity by SDS Page	Conforms
> = 95%	
Reconstitution	Conforms
reconstitute in 250 ul sterile dH2O	

Specification: PRD.0.ZQ5.10000077702

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.