


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
MorDalPhos Pd G4

Product Number: 900276

Formula: C₄₄H₅₇N₂O₄PPdS
Formula Weight: 847.39 g/mol

TEST	Specification
------	---------------

Appearance (Color)	White to Yellow to Brown
Appearance (Form)	Powder or Crystals
¹ H NMR Spectrum	Conforms to Structure
³¹ P NMR Spectrum	Conforms to Structure
Carbon Content	56.1 - 68.6 %
Nitrogen Content	
Specification: Report Only	
Sulfur Content	
Specification: Report Only	
ICP Major Analysis	Confirmed
Confirms Palladium Component	

Specification: PRD.0.ZQ5.10000073624

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.