

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tetrabutylammonium fluoride hydrate - 98%

Product Number: 241512
CAS Number: 22206-57-1
MDL: MFCD00149980
Formula: C₁₆H₃₆FN · xH₂O
Formula Weight: 261.46 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Yellow
Appearance (Form)	Powder or Crystal or Chunk(s)
Infrared spectrum	Conforms to Structure
Titration with HClO ₄ Dried Basis	97.5 - 102.5 %
Water (by Karl Fischer)	13.0 - 17.0 %

Specification: PRD.2.ZQ5.10000010917

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.