


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Hydroxybenzyl alcohol - 99%

Product Number: H20601
CAS Number: 620-24-6
MDL: MFCD00004643
Formula: C₇H₈O₂
Formula Weight: 124.14 g/mol



TEST

Specification

Appearance (Color)	Conforms to Requirements
White to Off-White to Brown	
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder and/or Chunks	
Infrared spectrum	Conforms to Structure
GC (area %)	≥ 98.5 %
Solubility (Turbidity)	Clear to Slightly Hazy
50mg/ml, Ethanol	
Solubility (Color)	Colorless to Brown

Specification: PRD.0.ZQ5.10000011327

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.