

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

2-Oxo-1-imidazolidinecarbonyl chloride - 96%

Product Number:

524832

CAS Number:

13214-53-4

MDL:


MFCD02683529

Formula:

C₄H₅ClN₂O₂

Formula Weight:

148.55 g/mol

**TEST****Specification**

Appearance (Color)

White to Beige

Appearance (Form)

Powder or Crystals

Infrared spectrum

Conforms to Structure

Titration by AgNO₃

95.5 - 104.5 %

After Oxygen Combustion

Specification: PRD.0.ZQ5.10000029472

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.