

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Sodium hexafluoroantimonate(V) – technical grade

Product Number: 237981
CAS Number: 16925-25-0
MDL: MFCD00003482
Formula: F₆NaSb
Formula Weight: 258.74 g/mol



TEST

Specification

Appearance (Color)
Appearance (Form)
X-Ray Diffraction

White to Off-White
Powder or Chunks
Conforms to Structure

Specification: PRD.0.ZQ5.10000037976

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.