

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com

Email USA: techserv@sial.com


Outside USA: eurtechserv@sial.com

Product Specification

Product Name:
BES sodium salt - Vetec™ reagent grade

Product Number: **V900440**
CAS Number: 66992-27-6

Formula: C₆H₁₄NNaO₅S
Formula Weight: 235.2 g/mol



TEST	Specification
Appearance (Colour)	White
Appearance (Form)	Crystalline Powder
Solubility (Colour)	Colorless
Solubility (Turbidity)	Clear
333.3 mg/ml of H ₂ O	
Loss on Drying	≤ 1 %
Titration by H ₂ SO ₄ (Dry Basis)	≥ 99 %
Heavy Metals (as Lead)	< = 5 ppm

Specification: PRE.0.ZQ5.10000002552

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.