

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Glutathione-(glycine-¹³C₂,¹⁵N) trifluoroacetate salt - ≥99 atom % ¹³C, ≥98 atom % ¹⁵N, ≥95% (CP)

Product Number:

683620

MDL:

MFCD07369422



Formula Weight:

310.26 g/mol

Storage Temperature:

-20 °C

TEST	Specification
Appearance (Form)	Solid
Purity (HPLC)	≥ 95.0 %
¹³ C Enrichment	≥ 99 %
¹⁵ N Enrichment	≥ 98 %
LCMS	Conforms

Specification: PRD.0.ZQ5.10000101878

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.