


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Tetrahydrofuran solution - contains 20 % (v/v) 1,2-Propanediol

Product Number: **680303**



TEST	Specification
Appearance (Color)	Colorless
Appearance (Form)	Liquid
Color Test	≤ 10 APHA
Water (by Karl Fischer)	≤ 0.20 %
Titration Acid (meq/g)	≤ 0.02
Titration Base (meq/g)	≤ 0.02
Mass Balance	79 - 81 % (V)
V% Tetrahydrofuran	
Mass Balance	19 - 21 % (V)
V% 1,2-Propanediol	
Expiration Date Period	-----
2 years	

Specification: PRD.4.ZQ5.10000010612

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.