

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Monoclonal Anti-Myosin (Skeletal, Fast) antibody produced in mouse – clone MY-32, ascites fluid

Product Number:

M4276**TEST****Specification**

Appearance (Form)	Liquid
Appearance (Turbidity)	Clear to Slightly Hazy
Appearance (Color)	Very Faint Yellow to Very Dark Yellow to Very Dark Brown
Protein by Biuret	15 - 55 mg/ml
Mouse IgG1 by RID	5.0 - 10.0 mg/ml
WB-Cell Line/Tissue Extract	Rabbit Leg or Skeletal Muscle Myosin
WB Titer	1:1000-1:2000
IHC - Tissue	Rabbit Leg Muscle or Rabbit Tongue or Rat Tongue or Human Tonsil with Skeletal Muscle
IHC - Fixation	FFPE-Formalin Fixed Paraffin Embedded
IHC - Titer	1:400-1:800
Recommended Retest Period	-----
5 Years	

Specification: PRE.0.ZQ5.10000016212

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.