


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
(Bicyclo[2.2.1]hepta-2,5-diene)tetracarbonylmolybdenum(0) - 95%

Product Number: 119261
CAS Number: 12146-37-1
MDL: MFCD00065319
Formula: C₁₁H₈MoO₄
Formula Weight: 300.12 g/mol



TEST

Specification

Appearance (Color)

Conforms to Requirements

Yellow to Tan to Dark Green

Appearance (Form)

Powder

Proton NMR spectra

Conforms to Structure

Carbon

41.6 - 46.5 %

Specification: PRD.2.ZQ5.10000033584

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.