

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Aluminum – rod, diam. × L 3.0 mm × 100 mm, 99.999% trace metals basis

Product Number:

202576

CAS Number:

7429-90-5

MDL:

MFCD00134029

Formula:

Al

Formula Weight:

26.98 g/mol

Al

TEST	Specification
------	---------------

Appearance (Color)	Conforms
Silver	
Appearance (Form)	Conforms
Rods	
Diameter	Conforms
3 mm	
Length	Conforms
100 mm	
ICP Major Analysis	Confirmed
Confirms Aluminum Component	
Purity	Conforms
99.999% Based On Trace Metals Analysis	
Trace Metal Analysis	≤ 20.0 ppm

Specification: PRD.0.ZQ5.10000053369

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.