


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
5-Aminoindole - 97%

Product Number: **A59654**
CAS Number: 5192-03-0
MDL: MFCD00005679
Formula: C₈H₈N₂
Formula Weight: 132.16 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
Tan to Purple-Brown to Very Dark Brown	
Appearance (Form)	Conforms to Requirements
Powder, Crystals, Crystalline Powder, and/or Solid	
Proton NMR spectra	Conforms to Structure
Carbon	70.16 - 75.25 %
Nitrogen	20.46 - 21.94 %
Purity (Titration by HCL04)	96.5 - 103.5 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.