

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
GermananeProduct Number: **906026**
CAS Number: 13572-99-1**GeH**Formula: GeH
Formula Weight: 73,647 g/mol

TEST	Specification
Appearance (Color)	Grey to Black
Appearance (Form)	Powder
Infrared Spectrum	Conforms to Structure
X-Ray Diffraction	Conforms to Structure
Calcium (Ca)	≤ 1 %
Hydrogen Content	
Specification: Report Only	
Analysis Type	Conforms to Structure
X-ray photoelectron spectroscopy (XPS)	
Recommended Retest Period	-----
1 Year	

Specification: PRD.0.ZQ5.10000101455

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.