

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
2,1,3-Benzothiadiazole-4,7-bis(boronic acid pinacol ester) - 95%

Product Number: **702803**
CAS Number: 934365-16-9

Formula: C₁₈H₂₆B₂N₂O₄S
Formula Weight: 388.10 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White to Dark Yellow and Faint Orange to Dark Orange and Faint Red to Dark Red and Faint Beige to Dark Beige
Appearance (Form)	Conforms to Requirements
Solid, Powder, Crystals and/or Chunks	
Infrared spectrum	Conforms to Structure
Carbon	52.6 - 58.8 %
Nitrogen	6.8 - 7.6 %
Sulfur	7.8 - 8.7 %
Purity (GC)	≥ 94.5 %

Specification: PRD.1.ZQ5.10000024436

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.