


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
3-Thienylboronic acid - $\geq 95.0\%$

Product Number: **436844**
CAS Number: 6165-69-1
MDL: MFCD00151851
Formula: C₄H₅BO₂S
Formula Weight: 127.96 g/mol



TEST	Specification
Appearance (Color)	Conforms to Requirements
White to Tan	
Appearance (Form)	Conforms to Requirements
Crystals to Crystalline Powder and/or Chunks	
Proton NMR Spectrum	Conforms to Structure
Purity (Titration by NaOH)	$\geq 95.0\%$

Specification: PRD.0.ZQ5.10000038696

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.