


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Urapidil hydrochloride - solid

Product Number: U100
CAS Number: 64887-14-5
MDL: MFCD00078601
Formula: C₂₀H₂₉N₅O₃ · HCl
Formula Weight: 423.94 g/mol



TEST

Specification

Appearance (Color)	White to Off-White
Appearance (Form)	Powder
Solubility (Color)	Colorless to Very Faint Yellow
Solubility (Turbidity)	Clear
25 mg/ml, H ₂ O	
Infrared spectrum	Conforms to Structure
EmM	22 - 26
in MeOH	
Wavelength	267 - 269 nm
Purity (HPLC)	≥ 98 %

Specification: PRD.0.ZQ5.10000068446

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.