


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Neodymium(III) isopropoxide - 99.8% trace metals basis

Product Number: 410101
CAS Number: 19236-15-8
MDL: MFCD00058819
Formula: C₉H₂₁NdO₃
Formula Weight: 321.50 g/mol



TEST	Specification
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
ICP Major Analysis	Confirmed
Confirms Neodymium Component	
Purity	Meets Requirements
99.8% Based on Rare Earth Analysis	
Trace Rare Earth Analysis	≤ 2500.0 ppm

Specification: PRD.1.ZQ5.10000039939

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.