

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Phenoxyacetic acid sodium salt hydrate - ≥98%

Product Number:

194220

CAS Number:

313222-85-4

MDL:


MFCD00150795

Formula:

C₈H₇NaO₃ · xH₂O

Formula Weight:

174.13 g/mol

**TEST****Specification**

Appearance (Color)

White to Off-White

Appearance (Form)

Conforms to Requirements

Powder and/or Lumps

Infrared spectrum

Conforms to Structure

Purity (Titration by HClO₄)

98.0 - 102.0 %

Specification: PRD.0.ZQ5.10000020629

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.