

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)


## Product Specification

Product Name:

(E)-N-(2-Aminoethyl)-4-[2-[4-(3-azidopropoxy)phenyl]diazenyl]benzamide hydrochloride

Product Number:

771139



Formula:

C<sub>18</sub>H<sub>21</sub>N<sub>7</sub>O<sub>2</sub> · HCl

Formula Weight:

403.87 g/mol

**TEST****Specification**

Appearance (Color)

White to Very Dark Yellow

Appearance (Form)

Powder or Crystals

Proton NMR Spectrum

Conforms to Structure

Carbon

42.5 - 64.5 %

Nitrogen

19.3 - 29.3 %

Purity (HPLC)

≥ 79.5 %

Specification: PRD.1.ZQ5.10000035736

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.