

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Angiotensin Converting Enzyme from rabbit lung – ≥ 2.0 units/mg protein (modified Warburg-Christian)

Product Number:

A6778

CAS Number:

9015-82-1

MDL:

MFCD00130325

Storage Temperature:

-20 °C

TEST	Specification
------	---------------

Appearance (Color)

Off White to Yellow and Faint Beige to Beige

Appearance (Form)

Powder

units/mg protein

 ≥ 2.0

One unit will produce 1.0 micromole of hippuric acid from hippuryl-HIS-LEU per minute in 50 mM HEPES and 300 mM NaCl at pH 8.3 at 37 deg C.

% Protein (WC)

70 - 100

Specification: PRD.0.ZQ5.10000033241

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.