


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Polypropylene – Amorphous

Product Number: 428175
CAS Number: 9003-07-0
MDL: MFCD00084447
Formula: C₃H₆



TEST	Specification
Appearance (Color)	White to Off-White
Appearance (Form)	Conforms
Rubbery Pellets	
Infrared spectrum	Conforms to Structure
Softening Point	150.0 - 160.0 °C
Gardner Color	≤ 2
Viscosity	700 - 1400 cps
@ 190 degrees C	

Specification: PRD.2.ZQ5.10000046197

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.