

3050 Spruce Street, Saint Louis, MO 63103, USA  
Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)  
Email USA: [techserv@sial.com](mailto:techserv@sial.com)  
Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

Product Specification

Product Name:  
Anti-Thioredoxin antibody produced in rabbit – IgG fraction of antiserum, buffered aqueous solution

Product Number: T0803

Storage Temperature: -20 °C

| TEST                         | Specification          |
|------------------------------|------------------------|
| Appearance (Form)            | Liquid                 |
| Appearance (Turbidity)       | Clear to Slightly Hazy |
| Appearance (Color)           | Colorless              |
| IgG Concentration (OD 280nm) | 10.0 - 20.0 mg/ml      |
| DID-Titer                    | 1:8-1:16               |
| Dot Blot-Titer               | 1:5,000-1:10,000       |
| WB-Cell Line/Tissue Extract  | E. Coli                |
| WB Titer                     | 1:5,000-1:10,000       |
| Recommended Retest Period    | -----                  |
| 4 Years                      |                        |

Specification: PRE.0.ZQ5.10000019723

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.