

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:


Luperox® A75FP, Benzoyl peroxide, 75% remainder water – contains 25 wt. % water as stabilizer, 75%

Product Number:

513474

MDL:

MFCD00003071



Formula Weight:

242.23 g/mol

Storage Temperature:

2 - 8 °C

TEST	Specification
Appearance (Color)	White
Appearance (Form)	Powder
Infrared spectrum	Conforms to Structure
Assay	73.0 - 77.0 %
Active Oxygen	4.82 - 5.09 %
Vendor Information	Conforms
Product of Arkema, Inc.	
Recommended Retest Period	-----
1 Year	

Specification: PRD.1.ZQ5.10000019879

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.