

3050 Spruce Street, Saint Louis, MO 63103, USA


Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Anthraquinone – Vetec™ reagent grade

Product Number: **V900763**
CAS Number: 84-65-1

Formula: C₁₄H₈O₂
Formula Weight: 208.22 g/mol



TEST

Specification

Appearance (Colour)
Yellow to Light Green to Tan
Appearance (Form)
Infrared spectrum
Purity (HPLC)

Conforms to Requirements

Powder
Conforms to Structure
≥ 96.5 %

Specification: PRE.0.ZQ5.10000002515

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.