

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

11-Mercaptoundecyl trifluoroacetate - 99%

Product Number:

674230

CAS Number:

138524-05-7

MDL:

MFCD09842697

Formula:


C<sub>13</sub>H<sub>23</sub>F<sub>3</sub>O<sub>2</sub>S

Formula Weight:

300.38 g/mol

Storage Temperature:

2 - 8 °C



### TEST

### Specification

Appearance (Color)

Colorless to Yellow

Appearance (Form)

Liquid

Infrared spectrum

Conforms to Structure

Purity (HPLC)

≥ 98.5 %

Specification: PRD.0.ZQ5.10000054850

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.