

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaldrich.com](http://www.sigmaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)

## Product Specification

Product Name:

Poly(ethyleneimine) solution - average  $M_n$  ~60,000 by GPC, average  $M_w$  ~750,000 by LS, 50 wt. % in H<sub>2</sub>O

Product Number:


**181978**

CAS Number:

9002-98-6

MDL:

MFCD00084427

**TEST****Specification**

Appearance (Color)

Colorless to Light Yellow

Appearance (Form)

Viscous Liquid

Infrared spectrum

Conforms to Structure

Residue on Evaporation

48.0 - 52.0 %

Viscosity (mPa.s)

18000 - 40000

Specification: PRD.0.ZQ5.10000012160

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.