


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Bis(diphenylphosphino)acetylene - 98%

Product Number: 287954
CAS Number: 5112-95-8
MDL: MFCD00003045
Formula: C₂₆H₂₀P₂
Formula Weight: 394.38 g/mol



TEST	Specification
Appearance (Color)	White to Light Yellow
Appearance (Form)	Powder or Crystals
Infrared spectrum	Conforms to Structure
Carbon	77.2 - 81.1 %
Solubility (Turbidity) 25mg/ml, CH ₂ Cl ₂	Clear
Solubility (Color)	Colorless
Purity (TLC)	≥ 97 %

Specification: PRD.0.ZQ5.10000101816

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.