

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Gramicidin A from *Bacillus brevis*

Product Number: **50845**
CAS Number: 11029-61-1

Formula: C₁₄₈H₂₁₀N₃₀O₂₆
Formula Weight: 1.880 g/mol
Storage Temperature: 2 - 8 °C

TEST	Specification
------	---------------

Appearance (Form) Solid	Conforms
Appearance (Color) White to Off-White	Conforms
Purity (HPLC Area %)	≥ 90.0 %
HPLC Area minor component(% of area) Gramicidin C	≤ 5.0 %
Water Content	≤ 5.0 %

Specification: PRE.3.ZQ5.10000014230

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.