

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: [www.sigmaaldrich.com](http://www.sigmaaldrich.com)Email USA: [techserv@sial.com](mailto:techserv@sial.com)Outside USA: [eurtechserv@sial.com](mailto:eurtechserv@sial.com)**Product Specification**

Product Name:

2,4-Difluorobenzaldehyde 98%

Product Number:

265179

CAS Number:

1550-35-2

MDL:


MFCD00010326

Formula:

C<sub>7</sub>H<sub>4</sub>F<sub>2</sub>O

Formula Weight:

142.1 g/mol

**TEST****Specification**Appearance (Color)  
Appearance (Form)Colorless to Faint Yellow  
Liquid

Infrared spectrum

Conforms to Structure

Color Test

≤ 100 APHA

Purity (GC)

≥ 97.5 %

**Remarks:**

Specification Date : 12/16/2010

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at [Sigma-Aldrich.com](http://Sigma-Aldrich.com). For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.