


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
SPhos Pd G4

Product Number: 804282

Formula: C₄₀H₅₀NO₅PPdS
Formula Weight: 793.22 g/mol

TEST	Specification
Appearance (Color)	White to Yellow to Orange
Appearance (Form)	Powder or Crystals
¹ H NMR Spectrum	Conforms to Structure
³¹ P NMR Spectrum	Conforms to Structure
ICP Major Analysis	Confirmed
Confirms Palladium Component	
Carbon Content	57.1 - 63.8 %
Purity (HPLC)	≥ 95.5 %

Specification: PRD.2.ZQ5.10000077904

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.