

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

N-Boc-trans-4-hydroxy-L-proline methyl ester - 97%

Product Number:

639737

CAS Number:

74844-91-0

MDL:


MFCD00076981

Formula:

C₁₁H₁₉NO₅

Formula Weight:

245.27 g/mol

**TEST****Specification**

Appearance (Color)

White to Tan

Appearance (Form)

Conforms to Requirements

Powder, Crystals or Crystalline Powder and/or Chunks

Infrared spectrum

Conforms to Structure

Optical Rotation

-68.0 - -62.0 deg

c = 1%, CHCl₃

GC (area %)

≥ 96.5 %

Specification: PRD.0.ZQ5.10000014091

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.