


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Isopropenylboronic acid pinacol ester - contains phenothiazine as stabilizer, 95%

Product Number: 663212
CAS Number: 126726-62-3
MDL: MFCD08276843
Formula: C₉H₁₇BO₂
Formula Weight: 168.04 g/mol
Storage Temperature: 2 - 8 °C



TEST	Specification
------	---------------

Appearance (Color)	Colorless to Very Dark Green to Very Dark Brown
Appearance (Form)	Liquid
Infrared spectrum	Conforms to Structure
Purity (GC)	≥ 94.0 %
Excluding Stabilizer	

Specification: PRD.1.ZQ5.10000034677

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.