

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:

Strontium nitrate - 99.995% trace metals basis

Product Number:**204498**

CAS Number:

10042-76-9

MDL:

MFCD00011248

Formula:

N₂O₆Sr

Formula Weight:

211.63 g/mol



TEST	Specification
------	---------------

Appearance (Color)	White
Appearance (Form)	Crystals or Crystals with Chunks
Complexometric EDTA % Sr	40.8 - 42.0 %
ICP Major Analysis Confirms Strontium Component	Confirmed
Purity 99.995% Based On Trace Metals Analysis	Meets Requirements
Trace Metal Analysis	≤ 55.0 ppm

Specification: PRD.0.ZQ5.10000063799

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.