

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Ammonium carbonate - 99.999% trace metals basis

 Product Number:
 379999

 CAS Number:
 506-87-6

MDL: MFCD00010890

Formula: CH8N2O3
Formula Weight: 96.09 g/mol

HO OH 2NH

TEST

Appearance (Color)

Appearance (Form)

Crystals or Crystalline Powder

Infrared spectrum

Purity

99.999% Based On Trace Metals Analysis

Trace Metal Analysis

Specification: PRD.2.ZQ5.10000029732

Specification

White

Conforms to Requirements

Conforms to Structure Meets Requirements

≤ 15.0 ppm

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.